

2-[O-(2-Azidoethyl)-D-serine]-cyclosporin A

Chemical Properties

CAS No. :	164028-89-1
Formula:	C ₆₄ H ₁₁₄ N ₁₄ O ₁₃
Molecular Weight:	1287.7
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	2-[O-(2-Azidoethyl)-D-serine]-cyclosporin A (Compound 33) is a derivative of Cyclosporine. This compound exhibits an IC ₅₀ of 0.11 µg/mL in mixed lymphocyte reaction assays.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7766 mL	3.8829 mL	7.7658 mL
5 mM	0.1553 mL	0.7766 mL	1.5532 mL
10 mM	0.0777 mL	0.3883 mL	0.7766 mL
50 mM	0.0155 mL	0.0777 mL	0.1553 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481