

NPN12830

Chemical Properties

CAS No. : 22112-83-0
Formula: C52H38N4O8
Molecular Weight: 846.9
Storage: Keep away from direct sunlight
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	NPN12830 is a dye-based compound utilized in various photodynamic therapy research studies.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1808 mL	5.9039 mL	11.8078 mL
5 mM	0.2362 mL	1.1808 mL	2.3616 mL
10 mM	0.1181 mL	0.5904 mL	1.1808 mL
50 mM	0.0236 mL	0.1181 mL	0.2362 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481