

Bengamin hydrochloride

Chemical Properties

CAS No. :	1004-23-5
Formula:	C ₇ H ₁₀ ClNO
Molecular Weight:	159.61
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Bengamin (4-(Aminomethyl)phenol) hydrochloride is a pharmaceutical intermediate used in the synthesis of benzylguanidine derivatives, which exhibit antibacterial properties.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.2653 mL	31.3264 mL	62.6527 mL
5 mM	1.2531 mL	6.2653 mL	12.5305 mL
10 mM	0.6265 mL	3.1326 mL	6.2653 mL
50 mM	0.1253 mL	0.6265 mL	1.2531 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481