

## N-Dodecanoyl-L-homoserine lactone-3-hydrazone-fluorescein

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:   
Keep away from direct sunlight  
Store at -20°C  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	N-Dodecanoyl-L-homoserine lactone-3-hydrazone-fluorescein is a fluorescent labeled probe. It is utilized in detection studies of O-DDHSL and 3O-C12-HSL.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481