


CRBN ligand-774

Chemical Properties

CAS No. :	2888604-15-5
Formula:	C ₁₈ H ₂₃ ClN ₄ O ₂
Molecular Weight:	362.86
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	CRBN ligand-774 is an E3 ubiquitin ligase cereblon (CRBN) ligand used for recruiting the cereblon protein. It can be linked via a linker to a target protein ligand to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7559 mL	13.7794 mL	27.5588 mL
5 mM	0.5512 mL	2.7559 mL	5.5118 mL
10 mM	0.2756 mL	1.3779 mL	2.7559 mL
50 mM	0.0551 mL	0.2756 mL	0.5512 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481