

Tert-butyl n-methyl-N-[2-(piperidin-4-yl)ethyl]carbamate

Chemical Properties

CAS No. :	171049-32-4
Formula:	C13H26N2O2
Molecular Weight:	242.358
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Tert-butyl n-methyl-N-[2-(piperidin-4-yl)ethyl]carbamate is a PROTAC linker used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1261 mL	20.6305 mL	41.2609 mL
5 mM	0.8252 mL	4.1261 mL	8.2522 mL
10 mM	0.4126 mL	2.063 mL	4.1261 mL
50 mM	0.0825 mL	0.4126 mL	0.8252 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481