

CRBN ligand-692

Chemical Properties

CAS No. :	2925093-34-9
Formula:	C14H13N3O3
Molecular Weight:	271.28
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	CRBN ligand-692 is an E3 ubiquitin ligase cereblon (CRBN) ligand, utilized for the recruitment of the cereblon protein. This compound can be linked via a linker to a target protein ligand to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6862 mL	18.4311 mL	36.8623 mL
5 mM	0.7372 mL	3.6862 mL	7.3725 mL
10 mM	0.3686 mL	1.8431 mL	3.6862 mL
50 mM	0.0737 mL	0.3686 mL	0.7372 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481