

(1s,3s)-3-Amino-2,2,4,4-tetramethylcyclobutyl dihydrogen phosphate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	(1s,3s)-3-Amino-2,2,4,4-tetramethylcyclobutyl dihydrogen phosphate is a PROTAC linker, utilized for the synthesis of PROTAC molecules.
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481