

CRBN ligand-648

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	CRBN ligand-648 is a ligand for the E3 ubiquitin ligase cereblon (CRBN), utilized to recruit cereblon proteins. CRBN ligand-648 can be linked to a target protein ligand via a linker to form a PROTAC (proteolysis targeting chimera).
Targets(IC50)	Ligands for E3 Ligase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481