

12-Amino-1-dodecanol

Chemical Properties

CAS No. :	67107-87-3
Formula:	C ₁₂ H ₂₇ NO
Molecular Weight:	201.35
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	12-Amino-1-dodecanol is a PROTAC linker utilized in the synthesis of PROTAC molecules.
-------------	----------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9665 mL	24.8324 mL	49.6648 mL
5 mM	0.9933 mL	4.9665 mL	9.933 mL
10 mM	0.4966 mL	2.4832 mL	4.9665 mL
50 mM	0.0993 mL	0.4966 mL	0.9933 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481