

Boc-PEG4-I

Chemical Properties

CAS No. :	1291090-31-7
Formula:	C ₁₅ H ₂₉ O ₆
Molecular Weight:	432.3
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Boc-PEG4-I is a PROTAC linker utilized in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3132 mL	11.566 mL	23.1321 mL
5 mM	0.4626 mL	2.3132 mL	4.6264 mL
10 mM	0.2313 mL	1.1566 mL	2.3132 mL
50 mM	0.0463 mL	0.2313 mL	0.4626 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481