

## CRBN ligand-651

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	CRBN ligand-651 is an E3 ubiquitin ligase cereblon (CRBN) ligand utilized for recruiting cereblon proteins. This compound can be linked to a target protein ligand via a linker to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481