

Tos-PEG3-Bn

Chemical Properties

CAS No. : 84131-04-4
Formula: C₂₀H₂₆O₆S
Molecular Weight: 394.482
Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Tos-PEG3-Bn is a PROTAC linker utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.535 mL	12.6749 mL	25.3498 mL
5 mM	0.507 mL	2.535 mL	5.070 mL
10 mM	0.2535 mL	1.2675 mL	2.535 mL
50 mM	0.0507 mL	0.2535 mL	0.507 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481