

Boc-C-PEG3-I

Chemical Properties

CAS No. :	2162121-00-6
Formula:	C12H23IO5
Molecular Weight:	374.22
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Boc-C-PEG3-I is a PROTAC linker used in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6722 mL	13.3611 mL	26.7222 mL
5 mM	0.5344 mL	2.6722 mL	5.3444 mL
10 mM	0.2672 mL	1.3361 mL	2.6722 mL
50 mM	0.0534 mL	0.2672 mL	0.5344 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481