

## 14-Bromotetradecanoic acid

## Chemical Properties

CAS No. :	676-26-6
Formula:	C <sub>14</sub> H <sub>27</sub> BrO <sub>2</sub>
Molecular Weight:	307.27
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	14-Bromotetradecanoic acid serves as a PROTAC linker for synthesizing PROTAC molecules.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2545 mL	16.2723 mL	32.5447 mL
5 mM	0.6509 mL	3.2545 mL	6.5089 mL
10 mM	0.3254 mL	1.6272 mL	3.2545 mL
50 mM	0.0651 mL	0.3254 mL	0.6509 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481