

## CRBN ligand-859

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	CRBN ligand-859 is an E3 ubiquitin ligase cereblon (CRBN) ligand designed to recruit cereblon proteins. It can be linked to a target protein ligand using a linker to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481