

mmu-miR-139-3p inhibitor

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	mmu-miR-139-3p inhibitor is a complementary single strand fully modified with methoxy to the mature miRNA. This miRNA inhibitor binds specifically to mature miRNA, preventing it from complementary pairing with target genes and thereby inhibiting miRNA function. It is useful for loss-of-function studies of miRNA.
Targets(IC50)	MicroRNA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481