

2-(2-(tert-Butoxycarbonyl)-2-azaspiro[3.5]nonan-7-yl)acetic acid

Chemical Properties

CAS No. :	2306261-01-6
Formula:	C ₁₅ H ₂₅ NO ₄
Molecular Weight:	283.37
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	2-(2-(tert-Butoxycarbonyl)-2-azaspiro[3.5]nonan-7-yl)acetic acid is a PROTAC linker used in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.529 mL	17.6448 mL	35.2896 mL
5 mM	0.7058 mL	3.529 mL	7.0579 mL
10 mM	0.3529 mL	1.7645 mL	3.529 mL
50 mM	0.0706 mL	0.3529 mL	0.7058 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481