

## CRBN ligand-612

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

|               |  |
|---------------|--|
| Description   | CRBN ligand-612 is an E3 ubiquitin ligase cereblon (CRBN) ligand designed to recruit cereblon proteins. It can be linked to a target protein ligand via a linker to form a PROTAC. |
| Targets(IC50) | Ligands for E3 Ligase  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481