

## Ethyl 7-oxoheptanoate

## Chemical Properties

CAS No. :	3990-05-4
Formula:	C <sub>9</sub> H <sub>16</sub> O <sub>3</sub>
Molecular Weight:	172.22
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Ethyl 7-oxoheptanoate is a PROTAC linker used in the synthesis of PROTAC molecules.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8065 mL	29.0326 mL	58.0653 mL
5 mM	1.1613 mL	5.8065 mL	11.6131 mL
10 mM	0.5807 mL	2.9033 mL	5.8065 mL
50 mM	0.1161 mL	0.5807 mL	1.1613 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481