

3-Aminobicyclo[1.1.1]pentane-1-carbonitrile

Chemical Properties

CAS No. :	1936643-30-9
Formula:	C ₆ H ₈ N ₂
Molecular Weight:	108.14
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	3-Aminobicyclo[1.1.1]pentane-1-carbonitrile is a PROTAC linker, utilized in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	9.2473 mL	46.2364 mL	92.4727 mL
5 mM	1.8495 mL	9.2473 mL	18.4945 mL
10 mM	0.9247 mL	4.6236 mL	9.2473 mL
50 mM	0.1849 mL	0.9247 mL	1.8495 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481