

Ethyl hept-6-ynoate

Chemical Properties

CAS No. :	105339-02-4
Formula:	C ₉ H ₁₄ O ₂
Molecular Weight:	154.21
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Ethyl hept-6-ynoate is a PROTAC linker suitable for synthesizing PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.4847 mL	32.4233 mL	64.8466 mL
5 mM	1.2969 mL	6.4847 mL	12.9693 mL
10 mM	0.6485 mL	3.2423 mL	6.4847 mL
50 mM	0.1297 mL	0.6485 mL	1.2969 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481