

## 4-(1-Boc-4-piperidyl)phenylboronic acid pinacol ester

## Chemical Properties

CAS No. :	956136-85-9
Formula:	C <sub>22</sub> H <sub>34</sub> BNO <sub>4</sub>
Molecular Weight:	387.321
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	4-(1-Boc-4-piperidyl)phenylboronic acid pinacol ester is a PROTAC linker used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5818 mL	12.9092 mL	25.8184 mL
5 mM	0.5164 mL	2.5818 mL	5.1637 mL
10 mM	0.2582 mL	1.2909 mL	2.5818 mL
50 mM	0.0516 mL	0.2582 mL	0.5164 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481