

Fibronectin binding protein peptide, D3 peptide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fibronectin binding protein peptide, D3 peptide, originates from the D motif of the Staphylococcus aureus fibronectin receptor. Fibronectin binding protein peptide, D3 peptide, can inhibit the binding of 125I-Fn to <i>S. dysgalactiae</i> .
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481