

DMTr-2'-O-C6-rG(iBu)-3'-CE-phosphoramidite

Chemical Properties

CAS No. :	2642321-69-3
Formula:	C50H66N7O9P
Molecular Weight:	940.09
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	DMTr-2'-O-C6-rG(iBu)-3'-CE-Phosphoramidite is a monomeric material utilized in the synthesis of nucleic acids.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0637 mL	5.3186 mL	10.6373 mL
5 mM	0.2127 mL	1.0637 mL	2.1275 mL
10 mM	0.1064 mL	0.5319 mL	1.0637 mL
50 mM	0.0213 mL	0.1064 mL	0.2127 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481