

4-(Azetidin-3-yl)piperazin-2-one hydrochloride

Chemical Properties

CAS No. :	1624260-62-3
Formula:	C7H14ClN3O
Molecular Weight:	191.66
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	4-(Azetidin-3-yl)piperazin-2-one hydrochloride is a PROTAC linker that serves in the synthesis of PROTAC molecules.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2176 mL	26.0879 mL	52.1757 mL
5 mM	1.0435 mL	5.2176 mL	10.4351 mL
10 mM	0.5218 mL	2.6088 mL	5.2176 mL
50 mM	0.1044 mL	0.5218 mL	1.0435 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481