

Desamino lenalidomide-piperidine-C-OH

Chemical Properties

CAS No. :	2229717-72-8
Formula:	C19H23N3O4
Molecular Weight:	357.41
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Desamino lenalidomide-piperidine-C-OH is an E3 ligase ligand-linker conjugate that includes a CRBN-based ligand and a linker, suitable for the synthesis of PROTAC.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7979 mL	13.9895 mL	27.9791 mL
5 mM	0.5596 mL	2.7979 mL	5.5958 mL
10 mM	0.2798 mL	1.399 mL	2.7979 mL
50 mM	0.056 mL	0.2798 mL	0.5596 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481