

Rhodamine B-Chitosan (MW 200000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Rhodamine B-Chitosan (MW 200000) is a chitosan marked with Rhodamine B, possessing a molecular weight of 200000. It combines the polysaccharide properties of chitosan with the fluorescent characteristics of Rhodamine B, playing a crucial role in fields such as fluorescence labeling, cell imaging, and drug delivery research (Ex/Em = 546/610 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481