

## 3-(tert-Butyldimethylsilyloxy)propan-1-amine

## Chemical Properties

CAS No. :	115306-75-7
Formula:	C <sub>9</sub> H <sub>23</sub> NOSi
Molecular Weight:	189.37
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	3-(tert-Butyldimethylsilyloxy)propan-1-amine is a PROTAC linker utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2807 mL	26.4033 mL	52.8067 mL
5 mM	1.0561 mL	5.2807 mL	10.5613 mL
10 mM	0.5281 mL	2.6403 mL	5.2807 mL
50 mM	0.1056 mL	0.5281 mL	1.0561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481