

Sp-2'-O-Me-CTPαS tetrasodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sp-2'-O-Me-CTPαS tetrasodium is a CTP analogue. It is utilized in enzymatic processes to incorporate thiophosphate and a modified 2'-OH group into RNA.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481