

Caloxetate trisodium

Chemical Properties

CAS No. :	207230-20-4
Formula:	C ₂₃ H ₂₈ CaN ₃ Na ₃ O ₁₁
Molecular Weight:	631.53
Storage:	Keep away from direct sunlight Store at -20°C <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Gadoxetate disodium is a liver-specific magnetic resonance imaging (MRI) contrast agent used to detect and differentiate liver tumors.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5835 mL	7.9173 mL	15.8346 mL
5 mM	0.3167 mL	1.5835 mL	3.1669 mL
10 mM	0.1583 mL	0.7917 mL	1.5835 mL
50 mM	0.0317 mL	0.1583 mL	0.3167 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481