

## Ethyl 3-(4-(hydroxymethyl)phenyl)propanoate

## Chemical Properties

CAS No. :	107859-98-3
Formula:	C <sub>12</sub> H <sub>16</sub> O <sub>3</sub>
Molecular Weight:	208.26
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Ethyl 3-(4-(hydroxymethyl)phenyl)propanoate is a PROTAC linker used in the synthesis of PROTAC molecules.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8017 mL	24.0085 mL	48.0169 mL
5 mM	0.9603 mL	4.8017 mL	9.6034 mL
10 mM	0.4802 mL	2.4008 mL	4.8017 mL
50 mM	0.096 mL	0.4802 mL	0.9603 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481