

## 2-tert-Butyl 6-methyl 2-azaspiro[3.3]heptane-2,6-dicarboxylate

## Chemical Properties

CAS No. :	1408074-81-6
Formula:	C <sub>13</sub> H <sub>21</sub> NO <sub>4</sub>
Molecular Weight:	255.31
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	2-tert-Butyl 6-methyl 2-azaspiro[3.3]heptane-2,6-dicarboxylate is a PROTAC linker used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9168 mL	19.584 mL	39.1681 mL
5 mM	0.7834 mL	3.9168 mL	7.8336 mL
10 mM	0.3917 mL	1.9584 mL	3.9168 mL
50 mM	0.0783 mL	0.3917 mL	0.7834 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481