

Methyl 7-aminoheptanoate hydrochloride

Chemical Properties

CAS No. :	17994-94-4
Formula:	C ₈ H ₁₈ ClNO ₂
Molecular Weight:	195.69
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Methyl 7-aminoheptanoate hydrochloride is a PROTAC linker utilized in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.1101 mL	25.5506 mL	51.1012 mL
5 mM	1.022 mL	5.1101 mL	10.2202 mL
10 mM	0.511 mL	2.5551 mL	5.1101 mL
50 mM	0.1022 mL	0.511 mL	1.022 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481