

## CRBN ligand-469

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	CRBN ligand-469 is a ligand for the E3 ubiquitin ligase cereblon (CRBN), intended to recruit cereblon proteins. It can be linked via a linker to a target protein ligand to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481