

## CY7-Alginate

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:   
Keep away from direct sunlight  
Store at -20°C  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	CY7-Alginate is a CY7-labeled sodium alginate, which merges the polysaccharide properties of sodium alginate with the fluorescent characteristics of CY7, making it significant in fields like fluorescent labeling, cell imaging, and drug delivery research (Ex/Em = 740/770 nm).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481