

## Lacosamide impurity 1

## Chemical Properties

CAS No. :	171623-03-3
Formula:	C <sub>12</sub> H <sub>16</sub> N <sub>2</sub> O <sub>3</sub>
Molecular Weight:	236.27
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Lacosamide impurity 1 (Compound 2 in Figure 2 on page 6) is an impurity of Lacosamide and also serves as a key intermediate in its synthesis.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2324 mL	21.1622 mL	42.3245 mL
5 mM	0.8465 mL	4.2324 mL	8.4649 mL
10 mM	0.4232 mL	2.1162 mL	4.2324 mL
50 mM	0.0846 mL	0.4232 mL	0.8465 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481