

Aldoview

Chemical Properties

CAS No. :	1356479-94-1
Formula:	C ₁₅ H ₁₀ F ₃ N ₃
Molecular Weight:	289.26
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	AldoView (Compound 3) is an intermediate that can be labeled with fluorine-18 to synthesize [18F]AldoView. This compound serves as a highly selective aldosterone synthase PET imaging agent. AldoView is utilized in imaging studies for primary aldosteronism.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4571 mL	17.2855 mL	34.571 mL
5 mM	0.6914 mL	3.4571 mL	6.9142 mL
10 mM	0.3457 mL	1.7285 mL	3.4571 mL
50 mM	0.0691 mL	0.3457 mL	0.6914 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481