

C3 BODIPY Ganglioside GM1 (d18:1/3:0) ammonium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | C3 BODIPY Ganglioside GM1 (d18:1/3:0) ammonium is a fluorescently labeled GM1 ganglioside probe, useful as a biological experimental marker. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481