

Peptidyl Arginine Deiminase, Rabbit

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Peptidyl Arginine Deiminase, Rabbit (EC 3.5.3.15) is an enzyme that catalyzes the conversion of arginine to citrulline. Calcium ions are necessary for the activity of Peptidyl Arginine Deiminase.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481