

Solvent Blue 35-d17

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Solvent Blue 35-d17 (Sudan Blue II-d17) is the deuterium-labeled version of Solvent Blue 35. Solvent Blue 35 (Sudan Blue II; Oil Blue 35) is a blue solvent dye utilized for marking specific diesel fuels, among other applications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481