

## Malic Dehydrogenase, Porcine

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Malic Dehydrogenase, Porcine (EC 1.1.1.37) exists in eukaryotic cells as two isoenzymes: one is expressed in the mitochondria and plays a role in the tricarboxylic acid cycle, while the other is expressed in the cytoplasm, converting malate back to oxaloacetate from the mitochondria.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481