

Oxindole-d4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Oxindole-d4 (2-Indolinone-d4) is a deuterium-labeled version of Oxindole. Oxindole (Indolin-2-one) serves as an aromatic heterocyclic building block, with its derivatives emerging as key compounds in kinase inhibitor research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481