

L-Hyoscyamine-d3-1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	L-Hyoscyamine-d3-1 (Hyoscyamine-d3) is a deuterium-labeled variant of L-Hyoscyamine. L-Hyoscyamine (Daturine) is a naturally occurring tropane alkaloid found in plants and acts as an effective competitive antagonist of muscarinic receptors (MR). It is the levorotatory isomer of Atropine.
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481