

## Velnacrine-d3 maleate

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Velnacrine-d3-1 maleate (HP 029-d3-1) is the deuterium-labeled form of Velnacrine maleate. Velnacrine maleate (HP 029) is an orally active cholinesterase inhibitor utilized in Alzheimer's disease research.
Targets(IC50)	Cholinesterase (ChE)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481