

6-Phosphogluconic Dehydrogenase, Sheep

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	6-Phosphogluconic Dehydrogenase, Sheep (EC 1.1.1.44) is an enzyme in the oxidoreductase family that acts on the CH-OH group of donor molecules using NAD ⁺ or NADP ⁺ as an acceptor. The enzyme's substrates are 6-phospho-D-gluconate and NADP ⁺ , and its products are D-ribulose-5-phosphate, CO ₂ , and NADPH.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481