

(Rac)- $\delta$ -Viniferin

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

(Rac)- $\delta$ -Viniferin is the racemic form of  $\delta$ -Viniferin and an isomer of  $\epsilon$ -Viniferin (epsilon-Viniferin).  $\epsilon$ -Viniferin is a resveratrol dimer that exhibits potent inhibitory effects on all CYP enzyme activities with  $K_i$  values ranging from 0.5 to 20  $\mu$ M. It also offers significant antioxidant, anti-inflammatory, antidiabetic, and neuroprotective benefits.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481