

TEAD ligand-Linker Conjugate 1

Chemical Properties

CAS No. :	2918762-34-0
Formula:	C31H31ClF2N2O5
Molecular Weight:	585.04
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	TEADligand-Linker Conjugate 1 is a synthetic target protein ligand-linker conjugate, designed for the synthesis of PROTACs like PROTACTEADdegrader-2. PROTACTEADdegrader-2 serves as a potent TEAD1 PROTAC degrader with anti-cancer properties.
Targets(IC50)	Target Protein Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7093 mL	8.5464 mL	17.0928 mL
5 mM	0.3419 mL	1.7093 mL	3.4186 mL
10 mM	0.1709 mL	0.8546 mL	1.7093 mL
50 mM	0.0342 mL	0.1709 mL	0.3419 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481