

## (R)-Ureidoisobutyric acid

## Chemical Properties

CAS No. :	1315050-26-0
Formula:	C <sub>5</sub> H <sub>10</sub> N <sub>2</sub> O <sub>3</sub>
Molecular Weight:	146.15
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	(R)-Ureidoisobutyric acid is a metabolite of nucleosides.
Targets(IC50)	Nucleoside Antimetabolite/Analog

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.8423 mL	34.2114 mL	68.4229 mL
5 mM	1.3685 mL	6.8423 mL	13.6846 mL
10 mM	0.6842 mL	3.4211 mL	6.8423 mL
50 mM	0.1368 mL	0.6842 mL	1.3685 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481