

## EZH2 ligand-Linker Conjugate 1

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	EZH2ligand-Linker Conjugate 1 is a target protein ligand-linker conjugate that includes a ligand capable of binding to EZH2. It can be utilized for the synthesis of PROTAC EZH2 Degradar-9.
Targets(IC50)	Histone Methyltransferase,Ligands for Target Protein for PROTAC

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481